

Boardcare

Boardcare WiO reagent for the determination of water content in lubricants/marine liquid fuels



Contents:

Boardcare Water	20g x 4 tubes
Reagent Paste	
Boardcare WiO	1 x 500ml
Releaser	
250ml Plastic	1 Number
Bottle	
10ml Syringe	2 Number
Gloves	4 Pairs
50ml Graduated	1 Number
Tube	
Magnetic Pin,	1 Each
Steel Bal	

Order# H00002

One tube contains 20g of paste. When mixed with 125ml of Boardcare WiO releaser, it can be used for fourteen tests. One refill pack is sufficient for 50 tests if handled properly. For repeatable results, mix the paste with the liquid releaser just before testing. Avoid long-term storage of the ready-to-use liquid.

Boardcare water-in-oil test reagents can also be used with non-Boardcare brand test devices that follow the same operation technique.